


MEDIVATORS™
VALVE BASKET™
Single Use Cage



EASY, SAFE AND EFFECTIVE REPROCESSING OF YOUR REUSABLE ENDOSCOPE VALVES

With its simple twist-loop design, the Single-Use Valve Basket attaches your reusable valves to the endoscope and keeps each unique set of equipment together throughout the reprocessing cycle. The Single-Use Valve Basket's patented design ensures safe and effective decontamination, representing a major step forward compared to "tea strainers" and metal baskets, which can become a haven for debris and bio burden.

- Inexpensive, lightweight and easy to store
- Enables effective flow of detergents, disinfectants and water around the valves
- Click-lock catch secures valves to prevent inadvertent re-use

Valves (including rinsing valves) should stay with a named endoscope as a set, to prevent cross infection and enable full traceability [23], pg. 951

Beilenhoff U et al. ESGE-ESGENA guideline: Cleaning and disinfection in gastrointestinal endoscopy... Endoscopy 2008; 40: 939-957

REPROCESS REUSABLE VALVES IN 3 EASY STEPS

1

Place valves into basket



2

Close valve basket



3

Attach basket to endoscope umbilical cord



Ordering Information

Product Description	Part Number	Units Per Box
Single-use Valve Basket for decontamination of reusable endoscope valves	105025	100

VALVE BASKET™ is a trademark of Medivators Inc.



Manufactured in the UK by:

Cantel Medical (UK) Limited
Campfield Road
Shoeburyness
Essex
SS3 9BX
UK
Tel: +44 (0) 1702 291878

Cantel Medical (UK) Limited
Wolseley House
Staffordshire Technology Park
Beaconsdale, Stafford ST18 0GA
United Kingdom
Tel: +44 (0) 1785.782420
Fax: +44 (0) 1785.782421

Medivators BV
Sourethweg 11
6422 PC Heerlen
The Netherlands
Tel: +31.45.5.471.471
Fax: +31.45.5.429.595

Medivators Inc.
14605 28th Avenue North
Minneapolis, MN
55447-4822 USA
Tel: +1.763.553.3300
Fax: +1.763.553.3387

Cantel Medical Asia/Pacific Pte. Ltd.
1A International Business Park
#05-01 Singapore 609933
Tel: +65.6227.9698
Fax: +65.6225.6848

Cantel Medical Devices (China) Co. Ltd.
Unit 804-805, Innov Tower Block A,
Hongmei Road, Xuhui 200233 Shanghai
Tel: +86 21 60161380
Fax: +86 21 61210913