

MEDIVATORS®
RapidAER®
Endoscope Reprocessor



Extremely fast
17 minute cycle

Single bay for
maximum flexibility

Advanced pass
through design

EVOLUTIONARY PROCESS, REVOLUTIONARY SOLUTION

For over 20 years, Cantel Medical's endoscope decontamination solutions have delivered unparalleled levels of efficiency, reliability and economy to over 300 endoscopy units worldwide.

The MEDIVATORS® RAPIDAER® Endoscope Reprocessor offers a totally flexible solution to suit all sizes of endoscopy department and uses the latest modular construction methods to produce the fastest and smallest pass-through washer disinfector in the world.





MULTI RAPIDAER® ENDOSCOPE REPROCESSOR INSTALLATION

- Totally flexible installation options to suit any specification
- Easy to add more units as facility capacity grows
- All service access is from the dirty side minimising installation footprint

THE FASTEST WAY TO INCREASE PRODUCTIVITY

- Complete disinfection cycle times from 17 minutes
- Unique modular design for uptime increase to 95%
- External scope loading for efficient throughput
- Small footprint for maximum room capacity



THE SIMPLEST WAY TO ENSURE RELIABILITY

- Single bay design reduces machine downtime risk
- 90% of components are proven standard parts
- 95% of components accessed from loading side
- Intuitive, easy to use controls reduce cycle failures



THE INTELLIGENT WAY TO GUARANTEE SECURITY

- Automated daily Thermal Self Disinfection cycle at >80°C
- RO water used throughout each cycle
- Complete traceability through latest RFID technology
- Fully compliant with EN-15883



Product Specifications RAPIDAER® Endoscope Reprocessor



Specification

Model	RAPIDAER Endoscope Reprocessor
Dimensions	600mm (w) x 800mm (d) x 1930mm (h)
Processing Time	Approx. 17 minutes, dependant on water temperature
Disinfectant	RAPISPOR® two part, single shot, Peracetic acid based disinfectant
Detergent	INTERCEPT® PLUS alkaline single shot detergent
Power Requirements	230v, 50hz, 32 amp supply
Connection to supply	IEC 60309 industrial coupler
Construction	Steel frame, polycarbonate doors, plastic cover panels
Noise Level	< 58 dBA
Weight	240 kg
Process Capacity	1 large flexible endoscope per cycle
Operating Temperature Range	Ambient (0-40°C)
Humidity	20%-90%

Consumables

Description	Part Numbers
Air Extract Carbon Filter	IS-LA 7014
Printer Rolls (Pack of 20)	IS 795732
INTERCEPT® PLUS Detergent	900540
RAPISPOR® (2 x 5 litre bottles)	ML02-0147
Autowipes Disinfectant Wipes (6 x 200 wipe drums)	900505



Cantel Medical (UK) Limited
 Wolsley House, Staffordshire Technology Park
 Beaconside, Stafford, Staffordshire, St18 0GA
 United Kingdom
 Tel: +44 (0) 1785 782420
 Fax: +44 (0) 1785 782421

MEDIVATORS Inc. Headquarters
 14605 28th Avenue North
 Minneapolis, MN
 55447-4822 USA
 Tel: +1 (0) 763 553 3300
 Fax: +1 (0) 763 553 3387

MEDIVATORS BV
 Sourethweg 11
 6422 PC Heerlan
 The Netherlands
 Tel: +31 (0) 45 5 471 471
 Fax: +31(0)45 5 429 595

MEDIVATORS Asia/Pacific Pte. Ltd
 1 A International Business Park
 #05-01 Singapore 609933
 Tel: +65 (0) 6227 9698
 Fax: +65 (0) 6225 6848

MEDIVATORS Beijing Representative Office
 Room 708, 7th Floor Kaiheng Center
 No.2 Chaoyangmennei Street
 Dongchang District, Beijing China 100010
 Tel: +8610(0) 6567 8446
 Fax: +8610(0) 6567 8445