

ENDO STRATUS™

Irrigation Pump and CO₂ Insufflator



Advanced CO₂ Insufflation for Faster Procedures and Better Patient Outcomes.



QUICKER PATIENT TURNAROUND*

- Improved cecal intubation rates¹
- Greater small bowel intubation depths²
- 38% reduction in nursing attention

HIGHER PATIENT SATISFACTION

- Less discomfort, pain and bloating^{3,4,1}
- Quicker recovery times⁴
- Decreased gas distension
- Absorbed 150 times faster than room air and is promptly eliminated via the lungs

Patient Pain and Discomfort CO₂ vs Room Air



THE COMPLETE CIRCLE OF PROTECTION

As the global vanguard in infection prevention, **only Cantel delivers the Complete Circle of Protection**, a full-value, proactive partnership dedicated to helping you remove risk, streamline operational efficiencies and optimise your success.

PROCEDURE Reducing the risk of patient cross contamination is at the forefront of infection prevention. Cantel innovates infection control products designed to improve patient outcomes, while increasing procedural efficiency.



ENDO STRATUS™ CO₂ INSUFFLATOR

Compatible with

- Tank (D or E) or DISS Wall CO₂
- All major GI Endoscopes

ENDO STRATUS™ IRRIGATION PUMP

Compatible with

- All major GI Endoscopes

Pump includes comfortable, universal foot pedal

Warm Water Irrigation Improves Your Visibility and Your Patient's Comfort

Adjustable water bottle holder with integrated water heater


Easy and quick to adjust flow rate

Compatible with ENDOGATOR™ Tubing and Connectors


Automatic prime button provides instant irrigation upon foot pedal depression



CO₂ CARTS



- Conveniently transport and store two CO₂ tanks.
- ENDO STRATUS™ Irrigation Pump and CO₂ Insufflation units stack neatly on top platform.
- Non-porous straps secure ENDO STRATUS Irrigation Pump and CO₂ Insufflation units for easy transportation between rooms.
- Basket provides additional storage.



ENDO STRATUS™ CARTS ORDERING INFORMATION

PART NUMBER DIMENSIONS

EGA-7030	ENDO STRATUS™ IV Style Cart
EGA-7032	ENDO STRATUS™ IV Style Cart with Dual Yoke
EGA-7031	ENDO STRATUS™ Box Style Cart
EGA-7033	ENDO STRATUS™ Box Style Cart with Dual Yoke



ENDO STRATUS™ PUMPS ORDERING INFORMATION

PART NUMBER DIMENSIONS

EGA-500	ENDO STRATUS™ Irrigation Pump
EGA-501	ENDO STRATUS™ Insufflator Unit
EGA-7012	Tank Yoke
EGA-7029	Dual Tank Yoke
EGA-7020	3 Foot Wall Hose (Low Pressure)
EGA-7023	9 Foot Wall Hose (Low Pressure)
EGA-7011	3 Foot Tank Hose (High Pressure)
EGA-7024	9 Foot Tank Hose (High Pressure)

1. Singh, R., Neo, E., & Nordeen, N. (2012, July). Carbon dioxide insufflation during colonoscopy in deeply sedated patients. *World Journal of Gastroenterology*, 18(25), 3250-3253.
2. Bretthauer, MD, PhD, M. (2010). Turning science into clinical practice – the case of carbon dioxide insufflation. *Endoscopy* 2010, 42(1), 1104-1105. doi:http://dx.doi.org/ 10.1055/s-0030-1255973
3. Bretthauer, MD, PhD, M. (2007, September). Carbon dioxide insufflation improves intubation depth in double-balloon enteroscopy: a randomized, controlled, double-blind trial. *Endoscopy* 2007, 39(10), 1064-1067. doi:DOI 0.1055/s-2007-966990
4. Lynch, MBA, BSN, RN, Hayes BSN, RN, CGRN, Buffum DNSc, APRN, CS. (2012). Insufflation Using CO₂ vs Room Air During Colonoscopy: Comparison of Patient Comfort, Recovery Time, and Nursing Resources [PowerPoint slides]. Veterans Affairs Medical Center San Francisco.

ENDO STRATUS™, ENDO SMARTCAP™ and ENDOGATOR™ are trademarks of Medivators Inc.

MEDIVATORS IS NOW CANTEL

www.cantelmedical.eu

Medivators BV | p: +31 45 5 471 471 (customer service)

© 2019 Cantel Medical Corp. P/N 50098-583D

