


CANEXIS™

Endoscopy Solution



CANEXIS™ Integrated Workflow Solution links data and decision-making to support the efficient operation of your critical hospital departments. The fully integrated functionality brings together differing departmental software into one powerful tool and offers seamless integration with your hospital information system.



A POWERFUL TOOL TO MANAGE AND ELEVATE ALL ASPECTS OF YOUR ENDOSCOPY DEPARTMENT


CANEXIS™ Endoscopy Solution was developed in close collaboration with endoscopy experts to create a high quality, user-friendly tool that gives you peace of mind that your facility is always executing best practices.

The cutting-edge software is vendor-neutral, flexible and configurable to support your workflow needs. Unparalleled quality management tools and the ability to achieve full traceability of your procedures and endoscopes makes CANEXIS Endoscopy Solution an invaluable tool to manage your endoscopy department.

CANEXIS™ Endoscopy Solution is a cutting-edge, user-friendly tool that supports the workflow in your endoscopy department.

PLAN. TRACK. ELEVATE.

CANEXIS™ Endoscopy Solution gives you visibility across your department in one powerful tool. The software allows you to monitor every step of your patient's experience in your department to support quality patient care. The traceability, ownership and accountability that the software delivers helps support patient safety and infection prevention.



PLANNING

Plan and schedule procedures according to your department requirements. Only a minimal set of data is required; the patient (which is automatically selected using your HIS), date and time, examiner and procedure type. The software allows you to reserve endoscopes for specific procedures and simplifies bed management enabling your department to run more efficiently.

PROCEDURE

In the procedure room you are able to easily register all data that is related to the procedure. View patient data and register each step of the patient journey, as well as your staff, endoscopes, CSSD sets and consumables including the bedside clean.


REPROCESSING

At a glance you are able to determine the real-time status and location of your entire endoscope inventory. Identify which endoscopes need to be reprocessed and where in your facility they are located. Track which endoscopes are out for repair and see trends over time to help you make repair vs. replace decisions. Not only that, but you will also have a record that every step was followed during reprocessing.



INFORMATION AT YOUR FINGERTIPS

The reporting tool available within CANEXIS™ Endoscopy Solution provides powerful analytics to help you make informed, evidence-based decisions about your endoscopy department. Reports can be easily generated for audits, accreditations, or infection-related recalls and can be created specific to your needs.



SERVICE & MAINTENANCE SUPPORT


Buying software requires careful consideration of the service and maintenance support that protects your investment. Cantel is known for outstanding service and maintenance.

With a Cantel service contract you have the peace of mind that your software is being maintained to the highest standard. A Cantel service contract for CANEXIS Integrated Workflow Solution will ensure the ongoing performance of your software and offers you the support you need including a dedicated customer support line, user/super-user and administrator training, software updates and upgrades.

We offer a comprehensive range of service and maintenance contracts to cover your needs and we will work with you to determine which plan is best for your facility. No matter which plan you choose, our team will support you to ensure the smooth running of your department.

Please contact us to find out more about the different levels of cover we provide.





THE COMPLETE CIRCLE OF PROTECTION

For more than 40 years Cantel has been delivering innovation and excellence in endoscope reprocessing. We offer a comprehensive line of endoscope reprocessing products including market-leading detergents and disinfectants, technologically advanced automated endoscope reprocessors, transportation systems, drying and storage cabinets and equipment tracking systems.

Our ecosystem of product and service offerings comes together in our **Complete Circle of Protection**, an infection prevention program designed to help you streamline reprocessing workflow, improve department efficiency and reduce the risk of infection.

CANEXIS™ is a trademark of Cantel Medical Corp., its affiliates or related companies.

www.cantelmedical.co.uk

TO PLACE AN ORDER

p: +44 (0)1702 291878 | e: info@cantelmedical.co.uk

© 2019 Cantel Medical Corp. ML-CX-001-UK-02

