

# ISACLEAN™

Detergent

# ISASPOR™

Single Shot Disinfectant/Sterilant



## ISACLEAN™ DETERGENT AND ISASPOR™ SINGLE SHOT DISINFECTANT/STERILANT— THE RIGHT FORMULA FOR THE BEST PERFORMANCE AND RESULT

As part of the overall solution for the reprocessing of endoscopes, Cantel Medical offers different solutions that have been specifically designed to be used with the MEDIVATORS™ ISA™ Endoscope Reprocessor, in order to satisfy all the needs of the users.

The cold chemical sterilant/high level disinfectant ISASPOR Single Shot Disinfectant/Sterilant and the ISACLEAN Detergent have been developed to guarantee a safe and reliable process together with high compatibility with all endoscope brands.

# ISACLEAN™

Detergent



- Active and tested against microbial biofilm
- 4 enzymes active against proteins, lipids/fats, carbohydrates, cellular components
- Unique formula thanks to the presence of ISAZONE™ Activator (patented new microbiological agent)
- Compatible with all brands of endoscopes thanks to neutral pH
- Specifically developed for manual cleaning or use in MEDIVATORS™ ISA™ Endoscope Reprocessor

# ISASPOR™

Single Shot Disinfectant/Sterilant



- Synergic action of ISAZONE™ Activator and Peracetic Acid (PAA)
- Sterilization claim according to UNI EN ISO 14937 at room temperature
- Low starting PAA and pH 6 guarantee high compatibility with medical devices
- Disinfection contact time: 3 minutes
- Single-use chemistry; aldehyde-free
- Clinically, microbiologically and toxicologically tested
- Specifically designed for use in the MEDIVATORS™ ISA™ Endoscope Reprocessor

## Ordering Information

Part Number	Description
500189	ISACLEAN™ Detergent (2 x 5L)
500187	ISASPOR™ Single Shot Disinfectant/Sterilant (10L Sol. A + B)
500188	ISASPOR™ Single Shot Disinfectant/Sterilant (5L Sol. A + B)

ISA™, ISASPOR™, ISAZONE™ and ISACLEAN™ are rtrademark of Cantel Medical (Italy) S.r.l.  
MEDIVATORS™ is a trademark of Medivators Inc.



 Cantel Medical (Italy) S.r.l.  
Via Laurentina, 169  
00071 Pomezia (RM) Italia  
Tel.: +39 06 9145399  
Fax: +39 06 9146099  
www.cantelmedical.it

Medivators Inc.  
14605 28th Avenue North  
Minneapolis, MN  
55447-4822 USA  
Tel: +1.763.553.3300  
Fax: +1.763.553.3387

Medivators BV  
Sourethweg 11  
6422 PC Heerlen  
The Netherlands  
Tel: +31.45.5.471.471  
Fax: +31.45.5.429.595

Cantel (UK) Limited  
Campfield Road  
Shoeburyness, Essex  
SS3 9BX,  
United Kingdom  
Tel: +44 (0) 1702.291878  
Fax: +44 (0) 1702.290013

Cantel Medical Asia/Pacific Pte. Ltd.  
1A International Business Park  
#05-01 Singapore 609933  
Tel: +65.6227.9698  
Fax: +65.6225.6848

Cantel Medical Devices (China) Co. Ltd.  
Unit 804-805, Innov Tower Block A,  
Hongmei Road, Xuhui 200233 Shanghai  
Tel: +86 21 60161380  
Fax: +86 21 61210913