

SCOPE BUDDY™

Endoscope Flushing Aid



CONSISTENT, HANDS-FREE FLUSHING OF ENDOSCOPES

REDUCES THE RISK OF REPETITIVE MOTION INJURIES

- Replaces repetitive and labour-intensive syringing of channels
- Simple, hands-free operation reduces muscle fatigue
- Provides constant fluid delivery through all endoscope channels

PROVIDES CONSISTENT FLUSHING

- Continuous, timed fluid delivery ensures consistent, repeatable endoscope flushing
- Audible and visual alerts upon completion
- Compatible with all major endoscope brands
- Works with your existing sink; compact size for countertop or wall mounting
- Meets or exceeds endoscope manufacturers requirements for channel flushing volumes

SHORTENS FLUSHING TIME

- Thoroughly flushes most endoscopes in 35 seconds
- Flushes multiple endoscope channels simultaneously
- Quickly and quietly delivers fluids through endoscope channels

INTERCEPT™ DETERGENT PROVIDES SUPERIOR ENDOSCOPE CLEANING

- Fast, one-minute contact time
- Low foaming, neutral pH formulation safe for endoscope use
- Cost-effective concentrated formula reduces per-use costs



THE COMPLETE CIRCLE OF PROTECTION

As the global vanguard in infection prevention, **only Cantel delivers the Complete Circle of Protection**, a full-value, proactive partnership dedicated to helping you remove risk, streamline operational efficiencies and optimise your success.



Visual display of endoscope flushing time

User-selectable time setting

Audible and visual alerts upon completion

One-touch start/stop operation

Simple and quick tubing connections

*SCOPE BUDDY is only intended to be used as a Flushing Aid during Manual Cleaning step of Endoscope Reprocessing. It is not intended to manually clean endoscopes.

MANUAL CLEAN Manual cleaning is essential in removing bioburden and preparing the endoscope for high-level disinfection. For optimal results, Cantel's manual cleaning products and detergents are designed to be a critical first step to efficiently and effectively eliminating bacteria on endoscopes prior to reprocessing.

Complete endoscope connection kit



Specifications

Unit dimensions	24.1 H x 15.2 W x 20.3 D (cm) (9.5 H x 6 W x 8 D (inches))
Unit weight	4.5 kg (10 lbs)
Electrical requirements	110 – 230 V, 50/60 Hz
Power requirements	230v 50/60Hz
Electrical safety certifications	UL Standard 61010-1 CANADA/USA C22.2 No. 61010-1 ETL and cETL listed
Fluid delivery rate	1 litre/minute
Flush time	0 – 99 minutes
Timer control	Digital display with touchpad controls
System pressure	Built-in safety bypass
Cycle completion indicators	Audible and visual alerts
Unit placement	Counter-top or wall-mounting (bracket included)
Warranty	One year



Scope Buddy™ Flushing Aid information

PART NUMBER	DESCRIPTION
EFA-IN-G	SCOPE BUDDY™ Flushing Aid

Accessory ordering information

PART NUMBER	DESCRIPTION
78399-669	Complete endoscope connection kit

Detergent ordering information

PART NUMBER	DESCRIPTION
ML02-0134	INTERCEPT™ Detergent (4 x 3.8 litre)

SCOPE BUDDY™ and INTERCEPT™ are registered trademarks of Medivators Inc.

www.cantelmedical.co.uk

TO PLACE AN ORDER

p: 01702 291878 | e: orders@cantelmedical.co.uk

© 2019 Cantel Medical Corp. ML-SBD-001-UK-02

