

# TRANSCOPE & STORASCOPE™

Storage & Endoscope Transport Systems



# SAFE, SECURE TRANSPORT & STORAGE SYSTEMS FOR YOUR ENDOSCOPES

With their easy-to-manoeuvre and highly secure design, the TRANSCOPE™ & STORASCOPE™ Endoscope Storage and Transportation Systems protect your valuable equipment during transport to and from the procedure room, and while inventory remains in storage.

SURESTORE™  
Transportation Storage System



RAPIDAER™  
Endoscope Reprocessor



CLEANSCOPE™  
Transport & Short Term  
Storage System Tray



Conditioned Packaged  
Endoscope



STORASCOPE™ Cart



TRANSCOPE™ Cart

TRANSCOPE™ & STORASCOPE™  
Storage & Endoscope Transport Systems

“A lockable mobile trolley is advisable for the transportation of multiple endoscopes.”

BSG Guidance on Decontamination of Equipment for Gastrointestinal Endoscopy: 2017 Edition



## THE COMPLETE CIRCLE OF PROTECTION

As the global vanguard in infection prevention, **only Cantel delivers the Complete Circle of Protection**, a full-value, proactive partnership dedicated to helping you remove risk, streamline operational efficiencies and optimise your success.

## STORAGE AND ENDOSCOPE TRANSPORT SYSTEMS



The TRANSCOPE™ Cart accommodates 6 trays and is designed to safely and securely transport your clean or contaminated endoscopes and other delicate equipment to and from the decontamination facility.

The STORASCOPE™ Cart contains 10 trays and is designed for secure storage of valuable inventory within the decontamination/procedure room.

- 5" lockable castors for easy manoeuvrability and positioning of carts
- Locking doors for secure transport and storage
- Removable trays for easy cleaning

**DRY & STORE** Bacteria pose significant risk to endoscopes during transport and storage. Cantel's transport, drying and storage solutions are designed to protect valuable inventory, reduce cross contamination touchpoints, eliminate moisture in the endoscope channels and control humidity. Humidity or moisture in endoscope channels is known to aid bacterial growth.



## TRANSSCOPE™ & STORASCOPE™ Cart ordering information

PART NUMBER	DESCRIPTION	UNITS PER BOX
103543	TRANSSCOPE Cart incl. 6 Standard Trays	1
103544	TRANSSCOPE Cart incl. 6 Standard Trays & Tray Lids	1
103551	TRANSSCOPE Cart incl. 6 Contoured Trays	1
103555	TRANSSCOPE Cart incl. Zero Trays	1
103547	STORASCOPE Cart incl. 10 Standard Trays	1
103548	STORASCOPE Cart incl. 10 Standard Trays & Tray Lids	1
103553	STORASCOPE Cart incl. 10 Contoured Trays	1
103556	STORASCOPE Cart incl. Zero Trays	1

TRANSSCOPE™, STORASCOPE™, SURESTORE™ and CLEANASCOPE™ are trademarks of Cantel (UK) Limited

[www.cantelmedical.co.uk](http://www.cantelmedical.co.uk)

TO PLACE AN ORDER

p: 01702 291878 | e: [orders@cantelmedical.co.uk](mailto:orders@cantelmedical.co.uk)

© 2019 Cantel Medical Corp. ML-TS-001-UK-05

