

NEOCARE

Manual Decontamination Sinks | **Manual Clean**



in partnership with





THE COMPLETE CIRCLE OF PROTECTION

As the global vanguard in infection prevention, only Cantel delivers the Complete Circle of Protection, a full-value, proactive partnership dedicated to helping you remove risk, streamline operational efficiencies and optimise your success.

All Neocare products are constructed in the UK using EU certified 316-grade stainless steel and can be made to the size and layout that you require.

Key features:

- Foot controlled operation
- Free standing
- HTM 04 01 compliant
- Designed alongside NHS ergonomics experts to reduce stress on:
 - Wrists
 - Elbows
 - Shoulders
 - Neck
 - Back
- 300mm rear splashback to reduce splash contamination
- Anti-glare, easy clean finish
- Storage shelf beneath



Standard Pegler Mixer Taps

As standard the Manual Decontamination Sinks include:

- Two bowls and two work surfaces (or state your own configuration)
- Vertical height adjustment of 350mm (Standard operator height range circa 147cm to 199cm)
- Antibacterial feet
- Fully sealed waste outlet
- High-level weir overflow standpipe
- Pegler mixer taps
- Splashback located organisers to hold accessories per sink

Available at additional cost:

- Alternative taps, including HTM 64 and TMV3
- Deluge spray
- Jet spray gun
- Built-in chemical dosing
- Further organisers for storage


MANUAL CLEAN



Manual cleaning is essential in removing bioburden and preparing the endoscope for high-level disinfection. For optimal results, Cantel's manual cleaning products and detergents are designed to be a critical first step to efficiently and effectively eliminating bacteria on endoscopes prior to reprocessing.

NEOCARE

Manual Decontamination Sinks | **Manual Clean**



Product specifications Manual Decontamination Sink

| | |
|------------------------------------|--|
| Standard dimensions | 2000mm (w) x 650mm (d) x 850mm (starting height, vertical travel 350mm) (Can be manufactured to the size you require) |
| Standard sink bowl capacity | 32L |
| Construction | EU certified austenitic 316-grade stainless steel |
| Installation space | To allow free movement it is advisable to have a gap of 100mm to the left, right and rear of the unit. |

Manual Decontamination Sink ordering information

| Part number | Description | Units per box |
|-------------|-----------------------------|---------------|
| SINK | Manual Decontamination Sink | 1 |

Manual Decontamination Sink accessories ordering information

| Part number | Description | Units per box |
|-------------|--|---------------|
| TMTAP | Thermostatic Mixer Tap | 1 |
| TMTAPMB | Mono Block Thermostatic Mixer Tap | 1 |
| SELGUN1 | Selecta spray gun and fixings | 1 |
| SELGUN2 | Selecta spray gun, fixings and tools | 1 |
| DELUGESP | Deluge Spray | 1 |
| BRACKET1 | Brackets for Gloves, Buttons or Brushes | 1 |
| BRACKET2 | Brackets for Leak Testers, Flush and Chemical Dosing | 1 |
| BRACKET3 | Bracket for SCOPE BUDDY™ Endoscope Flushing Aid | 1 |

SCOPE BUDDY™ is a registered trademark of Medivators Inc.



Cantel (UK) Limited
Shoeburyness
Campfield Road
Essex, SS3 9BX
United Kingdom
Tel: +44 (0) 1702 291878
Fax: +44 (0) 1702 290013

Medivators BV
Sourethweg 11
6422 PC Heerlen
The Netherlands
Tel: +31.45.5.471.471
Fax: +31.45.5.429.595

Medivators Inc.
14605 28th Avenue North
Minneapolis, MN
55447-4822 USA
Tel: +1.763.553.3300
Fax: +1.763.553.3387

**TO PLACE
AN ORDER**

Call: 01702 291878
Fax: 01702 290013
Email: orders@cantelmedical.co.uk