

ENDO GATOR™

Irrigation Flushing Pump





THE COMPLETE CIRCLE OF PROTECTION

As the global vanguard in infection prevention, **only Cantel delivers the Complete Circle of Protection**, a full-value, proactive partnership dedicated to helping you remove risk, streamline operational efficiencies and optimise your success.

PROCEDURE Reducing the risk of patient cross-contamination is at the forefront of infection prevention. Cantel innovates infection control products designed to improve patient outcomes, while increasing procedural efficiency.

AUTOMATE IRRIGATION AND REDUCE PROCEDURE TIME

The ENDOGATOR™ Irrigation Flushing Pump eliminates manual syringing. Irrigation is provided via footswitch, reducing procedural time and increasing visibility.

SAFETY

- Reduced risk of cross-contamination
- Endoscope connectors also available in single patient version

COMPATIBILITY

- Works with all sterile water bottles
- Connects to all major brands of endoscopes

EFFICIENCY

- Minimises reprocessing requirements
- Creates consistent practices
- Cost-effective and completely disposable
- Endoscope connectors for all methods of irrigation (Biopsy, Jet)

Pump Specifications

Electrical

Voltage/Current Requirement:
120V/240V-50-60Hz

Pump










Flow Rate: 0-500 ml/min
Pressure: 0-48 psi

Ideal Flow Settings

Auxiliary Water Channel: 40-50
Biopsy Channel: 70 to Max

Physical Dimensions

Height: 152 mm Width: 178 mm
Depth: 203 mm Weight: 2.3 kg

	TUBING	BIOPSY CHANNEL IRRIGATION	JET CHANNEL IRRIGATION	
OLYMPUS® GI Endoscopes	 200230U	 100303	 100241	
PENTAX® GI Endoscopes		 100136	 100242	
STORZ / FUJIFILM™ GI Endoscopes	 200230FU	 100303	 Backflow Valve 100126	 Direct Luer Connection on Tubing

ENDOOGATOR™ and MEDIVATOR™ are trademarks of Medivators Inc.

All other company and product names are trademarks of their respective owner.

© 2017-2021 Cantel. All Rights Reserved. The Complete Circle of Protection diagram is a trademark owned by Cantel Medical Corp. and is subject to copyright protection owned by Cantel Medical Corp.

www.cantelmedical.eu

Medivators BV | t: +31 (0)45 547 1471 | e: csdistributors@medivators.com

© 2021 Cantel Medical Corp. LIT-00371/A-ENEU

